



PRP - quick, easy, safe


1.



CAUTION! First release plunger (A). Push towards tip!

- Draw blood into (A)
- Close device with white strile cap
- Invert 10 times, break off plunger

2.



CAUTION! Centrifuge break settings: slow or off

- Fill counterweight with 40 mL water
- Place (A) opposite counterweight in centrifuge
- Centrifuge 3 min. at 900 g (RCF)

3.



CAUTION! First release plunger (B). Push towards tip!

- Attach (B) to (A) with luer luer adapter
- Press (A) with the tapered end into the presser foot
- Transfer plasma into (B), stop at mark 1
- Attach air filter to (B) and draw back plunger completely

4.



- Close (B) with new white cap, break off plunger
- Align volume of counterweight to (B)
- Place (B) opposite counterweight in centrifuge
- Centrifuge 10 min. at 1300 g (RCF)

5.



- Attach 30 mL syringe to (B)
- Press (B) with the tapered end into the presser foot
- Extract PPP to 30 mL syringe to mark 4
- Exchange 30 mL syringe with 10 mL syringe

6.



- Transfer sterile air from 10 mL syringe to (B)
- Shake system for 1 minute and withdraw PRP
- Hold plunger of 10 mL syringe in place when detaching (B)
- Close 10 mL syringe until further use

