

The underlying principle of the HT Vista technology is that normal and malignant tissues display different heat transfer rates due to differences in composition, metabolism, tissue morphology, and vascular network, which affect their thermophysical properties. This innovative screening modality relies on the unique thermal signals recorded by the device, as the tissue is heated and left to cool down.

HT Vista is a proven technology tested on thousands of dogs. It's an easy-to-use, affordable point-of-care screening tool that will increase the number of diagnosed masses. A quick, non-invasive 40-second scan provides on-the-spot cloud analysis with results within 2-3 minutes!

98% The high Negative Predictive Value gives confidence in ruling out cancer

Promotion Code: VBSFREE13

Contact VBS Direct to arrange your In-Clinic Demonstration quoting the above code to secure your FIRST MONTH FREE (subject to terms and conditions please apply).



